**permobil** | SEATING + POSITIONING

**THE ADAPTOR<sup>®</sup> PAD  
and ROHO<sup>®</sup> DRY FLOATATION<sup>®</sup>  
Special Application Cushioning Devices**

Operation Manual

Products:

THE ADAPTOR® PAD (THE ADAPTOR)

ROHO® HEAL PAD® Cushion (HEAL PAD)

ROHO® MINI-MAX® Cushion (MINI-MAX)

ROHO® PACK-IT® Cushion (PACK-IT)

ROHO® Shower/Commode Cushion (Shower/Commode Cushion)

ROHO® Toilet Seat Cushion (Toilet Seat Cushion)

ROHO® Recliner System (Recliner System)

## Please Note

Supplier: This manual must be given to the user of this product.

Operator (Individual or Caregiver): Before using this product, read the instructions and save for future reference.

Products may be covered by one or more U.S. and foreign patents and trademarks, including ROHO®, DRY FLOATATION®, MINI-MAX®, PACK-IT®, HEAL PAD®, THE ADAPTOR®, Yellow Rope® and shape fitting technology®. Zytel® is a registered trademark of DuPont.

Product information is changed as needed; current product information is available at permobilus.com.

ROHO, Inc. is part of Permobil Seating and Positioning.

© 2007, 2019 Permobil

Rev. 2019-02-13

## Safety Symbols

### Warning:

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in serious injury.

### Caution:

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in minor to moderate injury or damage to the equipment or other property.

## Contact Information

**ROHO, Inc.**  
100 N. Florida Avenue  
Belleville, IL 62221-5429 U.S.A.

### U.S.A.:

800-851-3449  
Fax 888-551-3449  
Customer Support:  
orders.roho@permobil.com

### Outside the U.S.A.:

618-277-9150  
Fax 618-277-6518  
Customer Support:  
intlorders.roho@permobil.com

## Contents

Important Safety Information	2-3
Product Specifications	3
THE ADAPTOR HEAL PAD	4
MINI-MAX PACK-IT	5
Shower/Commode Cushion Toilet Seat Cushion	6
Recliner System	7
Cleaning and Disinfecting	8
Care Symbols Troubleshooting Storage, Disposal and Recycling Limited Warranty	9
Español	10
Français	19
Deutsch	28
Italiano	37
Nederlands	46
Dansk	55
Svenska	64
Norsk	73
Suomi	82
Português	91

## Important Safety Information

### Warnings:

- Check inflation frequently, at least once a day!
- Skin/soft tissue breakdown can occur due to a number of factors, which vary by individual. Check skin frequently, at least once a day. Redness, bruising, or darker areas (when compared to normal skin) may indicate superficial or deep tissue injury and should be addressed. If there is any discoloration to skin/soft tissue, STOP USE immediately. If the discoloration does not disappear within 30 minutes after disuse, immediately consult a clinician.
- The cushioning device and the cover MUST be compatible sizes and MUST be used as directed in this manual. IF THEY ARE NOT: 1) The benefits provided by the cushioning device may be reduced or eliminated, increasing the risk to skin or to other soft tissue, and 2) the individual may become unstable and vulnerable to falling.
- DO NOT use your cushioning device on top of, or in conjunction with, any other products or materials, except as indicated in this manual. Doing so may cause the individual to become unstable and vulnerable to falling.
- DO NOT place any obstructions between the individual and the cushion. USE ONLY a compatible cover with the cushion; placing any other item directly under the individual, whether inside or outside the cover, will make the cushion less effective.
- The cushion surface will acclimate to the surrounding temperature. Take appropriate precautions, especially when the cushion will contact unprotected skin.
- DO NOT use your cushioning device as a water flotation device (e.g., a Life Preserver). It will NOT support you in water.
- Use caution when skin will contact the seating surface. Tissue damage may result from increased friction.
- DO NOT use an under-inflated cushioning device or an over-inflated cushioning device because the product benefits will be reduced or eliminated, resulting in an increased risk to skin and other soft tissue. If the cushioning device does not appear to be holding air, or if you are not able to inflate or deflate the cushioning device, see "Troubleshooting". Immediately contact your clinical caregiver, distributor, or Customer Support if the problem persists.
- Use caution when handling the cushioning device: Metal grommets on the corners may be a finger trap hazard.
- Periodically, check the cover and components for damage, and replace if necessary.
- DO NOT expose the product to high heat, open flames or hot ashes. Testing or certification claims, including for flammability, may no longer apply to this device when it is combined with other products or materials. Check testing and certification claims for all products used in combination with this device.
- If you are unable to perform any tasks described in this manual, seek assistance by contacting your clinician, equipment provider, distributor, or Customer Support.

### Cautions:

- Keep your cushioning device away from sharp objects.
- DO NOT allow your cushioning device to come into contact with oil-based lotions or lanolin. They may degrade the material.
- Changes in altitude may require adjustment to a DRY FLOATATION cushioning device.
- Prolonged exposure to ozone may degrade materials used in the cushioning device, affect the performance of the product, and void the product warranty.

## Important Safety Information, continued

### Cautions, continued:

- DO NOT use an inflation valve as a handle for carrying or pulling your DRY FLOATATION cushioning device. Carry your cushioning device by the cushion base or the carrying handle on the cover.
- If a DRY FLOATATION cushioning device has been in temperatures lower than 32°F (0°C), and exhibits unusual stiffness, allow the cushioning device to warm to room temperature.
- Using a pump, cover or repair kit other than those compatible with your product may void your product's warranty.

*An experienced clinician should determine whether the product is appropriate for the individual's particular equipment needs. The product's intended use statement alone should not be relied upon to make that determination.*

*The medical products featured in this document are intended to be only one part of an overall care regimen that includes all therapeutic equipment and interventions. An experienced clinician should identify this care regimen after evaluating the individual's physical needs and overall medical condition.*

*A clinician should also evaluate for visual, reading, and cognitive impairment to determine the need for caregiver assistance or the need for other assistive technology, such as large-print instructions, to ensure proper product use.*

## Product Specifications

Package Contents: All products featured in this manual include an operation manual, product registration card, and Repair Kit. All products except THE ADAPTOR PAD include a ROHO Hand Pump. The MINI-MAX Cushion, PACK-IT Cushion, and Recliner Cushion include a cover.

To order replacement parts, contact your equipment provider or Customer Support.

Sizing Table	Height or Height Range		Width or Width Range		Depth or Depth Range		Average Weight*	
	in.	cm	in.	cm	in.	cm	lb.	kg
ADAPTOR	1	2	9	25	13½	34	½	0.2
HEAL PAD	½ - 1½	1.5 - 4	11½	29.5	12½	31.5	¾	0.3
MINI-MAX	1¼	3	15 to 20½	38.5 - 52	14¼ to 20¼	37.5 - 51.5	2¼	1
PACK-IT	1¼	3	16¼	41.5	9½	24	1½	0.7
Shower/Commode Seat	2¼	5.5	14¼ - 16¼	37.5 - 42	14¼ - 18¼	37.5 - 46.5	2	1
Toilet Seat	1 - 2¼	2.5 to 5.5	7¾	20	15¼	40	1¼	0.5
Recliner Cushion	Contact Customer Care for specifications.				* based on mid-sized cushioning devices			

Cushion expected lifetime: 5 years

The serial number label is located on the bottom of the cushion.

### Materials:

Cushioning devices: neoprene rubber (made without latex); Zytel valves

Components on some cushioning devices: Yellow Rope, buckles, snaps, hook and loop fastener, foam tab (on HEAL PAD)


Covers: Standard: two-way stretch top, spacer fabric sides. Heavy Duty: two-way stretch, fluid-resistant top and sides. Standard and Heavy Duty Covers: zipper or snaps and a non-skid bottom. Refer to the cover care label for material content and for compliance to flammability requirements. Elastic strap and metal ring included on the PACK-IT Cover.

**If provided, straps or hook and loop fasteners should be used to secure the cushioning device to the surface.**

## THE ADAPTOR

*Intended Use: THE ADAPTOR® PAD (THE ADAPTOR) is an air-filled, cellular-design, neoprene surface intended for cushioning and shock absorption needs and reduction of friction on surfaces such as a wheelchair side, headrest or footrest to provide skin/soft tissue protection. There is no weight limit. The ADAPTOR is not intended to be used as a supporting medium under the buttocks.*

**Instructions:** Use scissors to trim THE ADAPTOR to fit where needed (sides of wheelchairs, trays, head rests, foot rests). Cut between the sealed air cells when possible. To attach the product to the surface, use hook and loop fastener strips.


**Warning:** Air cells that have been cut **WILL NOT** provide cushioning.

## HEAL PAD

*Intended Use: The ROHO® HEAL PAD® Cushion (HEAL PAD) is an adjustable, air-filled, cellular-design, neoprene surface and is intended to conform to an individual's heel to provide skin/soft tissue protection and provide an environment to facilitate wound healing. There is no weight limit.*

### Inflation Instructions:


1. Place the HEAL PAD flat, with the cells facing up.
2. Slide the hand pump nozzle over the inflation valve. Open the inflation valve, at least one full rotation.
3. Inflate the HEAL PAD. Close the inflation valve until it comes to a stop. Remove the hand pump. Make sure the inflation valve is closed.


**Note:** The inflation valve will spin freely when open.

### Attachment Instructions:


- 1 Fasten the snaps.

- 2 Place the individual's heel on the cells.


- 3 Place the foam tab over the shin, insert the strap through the metal ring, and adjust the strap.

**CHECK INFLATION!** Confirm that there is **NO CONTACT** between the individual's heel and the base of the HEAL PAD. If there is contact, add air.

## MINI-MAX and PACK-IT

*Intended Use: The ROHO® MINI-MAX® Cushion (MINI-MAX) is an adjustable, air-filled, cellular-design, neoprene surface intended as a supplemental seating surface. There is no weight limit, yet the cushion must be properly sized to the individual.*


*Intended Use: The ROHO® PACK-IT® Cushion (PACK-IT) is an adjustable, air-filled, cellular-design, neoprene surface intended for short-term cushioning needs. There is no weight limit, yet the cushion must be properly sized to the individual.*

### Atmospheric Adjustment Instructions for MINI-MAX and PACK-IT:

1. Slide the hand pump nozzle over the inflation valve. Open the inflation valve, at least one full rotation. Inflate until all of the air cells feel firm.
2. Remove the hand pump. Do not close the inflation valve. Allow air to escape from the inflation valve.
3. Listen (or place your hand in front of the inflation valve) to confirm that deflation has stopped.
4. Close the inflation valve until it comes to a stop. Make sure the inflation valve is closed.
5. Place the cushion on the surface, making sure the air cells face the individual.

**Note:** *The inflation valve will spin freely when open.*

**CHECK INFLATION!** After atmospheric adjustment, confirm that there is **NO CONTACT** between the individual and the base of the cushioning device. If there is contact, add air.

When possible, use hook and loop fasteners to further secure the cushion to the seating surface. Hook fasteners are available upon request.

**Cover Removal:** Unzip the cover and remove the cushion.


**MINI-MAX Cover Replacement:** Insert the cushion into the cover. Pull the inflation valve and Yellow Rope through a corner opening.

**PACK-IT Cover Replacement:** Insert the PACK-IT into the cover. Close the zipper.

## Shower/Commode Cushion and Toilet Seat Cushion


*Intended Use: The ROHO® Shower/Commode Cushion (Shower/Commode Cushion) is an adjustable, air-filled, cellular-design, support surface that utilizes DRY FLOATATION® technology and is intended for short-term cushioning needs on bedside commodes and shower chairs to provide skin/soft tissue protection. There is no weight limit.*

**Shower/Commode Cushion**


*Intended Use: The ROHO® Toilet Seat Cushion (Toilet Seat Cushion) is an adjustable, air-filled, cellular-design, neoprene surface intended for short-term cushioning needs on toilets and portable bedside commodes. There is no weight limit.*

**Toilet Seat Cushion**

### Instructions for Shower/Commode Cushion and Toilet Seat Cushion:

1. Slide the hand pump nozzle over the inflation valve. Open the inflation valve, at least one full rotation. Inflate the cushion, and then close the inflation valve until it comes to a stop. Remove the hand pump. Make sure the inflation valve is closed. Repeat for all inflation valves.
2. Place the cushion on the surface, making sure the air cells face the individual. Confirm that the individual is in their usual position. Confirm correct placement of the cushioning device.
3. Adjust the inflation level as required. When adjustment is complete, close the inflation valve.

**Note: The inflation valve will spin freely when open.**

**Toilet Seat Cushion Attachment:** Use the straps and snaps to secure the cushion to the surface.

Contact between the skin and areas where the toilet seat is NOT covered by the Toilet Seat Cushion could result in tissue damage.

**CHECK INFLATION!** After adjustment, confirm that there is NO CONTACT between the individual and the base of the cushion. If there is contact, add air.

## Recliner System

*Intended Use: The ROHO® Recliner System (Recliner System) is an adjustable, air-filled, cellular-design, support surface that utilizes DRY FLOATATION® technology and is intended to conform to an individual's seated shape to provide localized skin/soft tissue protection and an environment to facilitate wound healing. There is no weight limit, yet the cushion must be properly sized to the individual.*

### Instructions:

1. Place the Recliner System on the chair with the air cells facing up and the non-skid surface of the cover facing the chair seat and chair back.


**When possible, use the hook or loop fasteners to further secure the cover to the seating surface. On the back cushion, hook and loop fasteners are located inside the cover and on the back of the cushion.**

2. Slide the hand pump nozzle over the inflation valve. Open the inflation valve, at least one full rotation. Inflate until all of the air cells feel firm. Close the inflation valve until it comes to a stop. Remove the hand pump. Make sure the inflation valve is closed. Repeat for all inflation valves.
3. Confirm that the individual is in their usual position. Confirm correct placement of the Recliner System.
4. Slide your hand between the cushion and the individual and feel for bony prominences. With your hand in place, open the inflation valve, releasing air as prescribed. Close the inflation valve. Repeat for all inflation valves.

**Note: The inflation valve will spin freely when open.**



**CHECK INFLATION!** Make sure that there is at least ½ in. (1.5 cm) of air between the cushion base and the individual's bony prominences. Follow the prescriber's recommended inflation levels.

Shown without cushion covers.

The seated figure shows standard orientation.

**Recliner Cushion Cover Removal:** Unsnap and remove the covers from the cushions.

### Recliner Cushion Cover Replacement:

1. The zipper on the Back Cushion cover should be closed. For the Seat Cushion and Back Cushion: The non-skid cover material should face up and the snapped opening of the cover should face you. The air cells should face down: Insert each cushion into its cover, inflation valve(s) first. Pull the inflation valve(s) and Yellow Rope through the corner opening(s). The inflation valves will protrude at the front of the Seat Cushion and at the top of the Back Cushion.
2. Place the two cushions with the snapped openings next to each other. Snap the Back Cushion to the Seat Cushion. Snap the bottoms of the covers together. Turn the cushions over so the air cells face up. Snap the tops of the covers together. The cushions should lie flat inside the covers.

## Cleaning and Disinfecting

Instructions apply to the products featured in this manual.

### Warnings:

- **Cleaning and disinfecting are separate processes. Cleaning must precede disinfection. Before use by a different individual: clean, disinfect, and check product for proper functioning.**
- **Follow all manufacturer bleach container and germicidal disinfectant safety guidelines.**
- **Note about sterilization: High temperatures accelerate aging and will damage the product assembly. The products featured in this manual are not packaged in a sterile condition, nor are they intended or required to be sterilized prior to use. If institution protocol demands sterilization: First, follow the cleaning and disinfecting instructions; then, open the inflation valve and use the lowest possible sterilizing temperature, but not to exceed 70°C (158°F), for the shortest time possible. DO NOT autoclave.**


**Cover Care:** Remove hook fastener pieces from bottom of cover prior to cleaning and disinfecting. For the cover on the Recliner System back cushion: Turn the cover inside out, close the zipper and snaps. Wash with like colors.

**To Clean the Cover:** First remove the cover from the cushioning device.

Machine wash in cold water (30°C) with mild detergent, gentle cycle, or wipe clean with neutral detergent and warm water. Tumble dry low.

**To Disinfect the Cover:** Hand wash in cold water (30°C) using 1 part household liquid bleach per 9 parts water, then rinse thoroughly. Or machine wash in hot water (60°C). Tumble dry low.

**Cushion Care Cautions:** DO NOT allow water or cleaning solution to enter the cushion. DO NOT machine wash or sanitize the cushion or machine dry the cushion. Use of the following may cause damage to the cushion: abrasives (steel wool, scouring pads); caustic, automatic dishwasher detergents; cleaning products that contain petroleum or organic solvents, including acetone, toluene, Methyl Ethyl Ketone (MEK), naphtha, dry cleaning fluid, adhesive remover; oil-based lotions, lanolin, or phenolic-based disinfectants; ozone gas; ultraviolet light cleaning methods.

**To Clean the Cushion:** Remove the cover and close the inflation valve. Hand wash, using mild liquid hand soap, hand dishwashing liquid, laundry detergent, or a multipurpose detergent (Follow product label instructions.) Gently scrub all cushion surfaces, using a soft plastic bristle brush, a sponge or a cloth. Rinse with clean water. Air dry. DO NOT expose to direct sunlight.

**Caution:** Rinse thoroughly. Detergent residue may cause air cells to stick together

**To Disinfect the Cushion:** Repeat cleaning instructions, using 1 part household liquid bleach per 9 parts water. Follow the safety guidelines on the bleach product label. Keep the product wet with the bleach solution for 10 minutes. Rinse thoroughly with clean water. Air dry. DO NOT expose to direct sunlight.

## Care Symbols

	Hand wash cold, maximum temperature shown.		Bleach (1 part bleach: 9 parts water)
	Machine wash hot, normal, maximum temperature shown.		Machine wash gentle, cold, maximum temperature shown.
	Hand wash.		Line dry in the shade.
			Tumble dry low.

## Troubleshooting

For additional assistance, contact your equipment provider, distributor, or Customer Support.

**Not Holding Air:** Inflate the cushion. Confirm that the inflation valve is completely closed. Look for holes in the cushion. If very small holes or no holes are visible, follow the instructions in the Repair Kit provided with your product. For inflation valve or hose damage, or for large holes or leaks in the cushion, see the Limited Warranty supplement, or contact Customer Support.

**Uncomfortable or Unstable:** Check the inflation level (see the product instructions). Confirm that the air cells and the top of the cover face the individual. Confirm that the cushioning device is not too big or too small for the surface. Allow the individual at least one hour to adjust to the cushioning device.

**Cushioning Device Slides:** Make sure the non-skid bottom of the cover is facing away from the individual. Confirm that the hook and loop fasteners, straps, or snaps are securely fastened. When possible, use hook and loop fasteners to further secure the cushion to the seating surface. Hook fasteners are available upon request.

## Storage, Disposal and Recycling

**Storage:** Clean and disinfect the product before storing. Open the inflation valve(s). Deflate the product (detach any product sections). Store the product in a container that protects the product from moisture, contaminants and damage.

**Disposal:** The components of the products in this manual are associated with no known environmental hazards when properly used and when disposed of in accordance with all local/regional regulations. Incineration must be performed by a qualified, licensed waste management facility.

**Recycling:** Contact a local recycling agency to determine recycling options for the product. For more recycling information, contact Customer Support.

## Limited Warranty

**Limited Warranty Term:** 24 months from the date the product was originally purchased. The warranty does not apply to punctures, tears, burns, or misuse, nor does it apply to the removable cover. See also the Limited Warranty supplement provided with your product, or contact Customer Support.

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MDSS GmbH